



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Samcomedical Hydraulic Stretcher Trolley 5 Functions

Product Code :

SAMCOMEDICAL10025A032

Description :

Samcomedical Hydraulic Stretcher Trolley 5 Functions

Technical Specification :

Samcomedical Hydraulic Stretcher Trolley 5 Functions with adjustable back rest, knee rest, trendelenburg, reverse trendelenburg & adjustable height. Back Rest adjustable by gas spring, knee rest adjustable by super smooth crank, other functions adjustable by hydraulic mechanism. Size: 2082 860 608/923 mm. Framework made of mild steel epoxy coated finish with X-ray permeable Top. Collapsible side railings. Equipped with hydraulic pump. Equipped with I.V. Pole, drainage hook & bumper protection. 150mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters



Samco Medical

Website: www.samcomedical.com,