



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Stretcher Trolley 4 Functions

Product Code :
SAMCOMEDICAL10025A033

Description :

Samcomedical Stretcher Trolley 4 Functions

Technical Specification :

Samcomedical Stretcher Trolley 4 Functions with adjustable back rest, trendelenburg, reverse trendelenburg & adjustable height. Back Rest adjustable by gas spring, other functions adjustable by hydraulic mechanism. Size: 1930 688 567/797 mm. Equipped with hydraulic pump. Equipped with I.V. Pole, drainage hook & bumper protection. 150mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. Two easy lifting guardrails (One on each side) made of PE engineering plastic, which are safe and reliable, and can be fixed upward and downward. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters



Samco Medical

Website: www.samcomedical.com,