



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Stretcher Trolley, 2 Functions

**Product Code :**  
SAMCOMEDICAL10025A035

**Description :**

Stretcher Trolley, 2 Functions

**Technical Specification :**

Stretcher Trolley, 2 Functions Stretcher trolley with adjustable Back rest & adjustable height. Back Rest adjustable by gas spring & Height adjustable by crank mechanism. Size: 1930  $\pm$  1/2 595  $\pm$  1/2 530/820 mm. Framework made of mild steel epoxy coated finish with ABS Plastic Top. Two easy lifting guardrails (One on each side) made of PE engineering plastic, which are safe and reliable, and can be fixed upward and downward. 125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. Equipped with separate directional wheel (5th Wheel) for easy operation only by one person. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),