



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Trolley

Product Code :
SAMCOMEDICAL10025A036

Description :

Stretcher Trolley

Technical Specification :

Stretcher Trolley Size: 1950L \pm 1/2 650W \pm 1/2 680 / 900 H mm. Height adjustable from 600 - 950 H mm. Frame work made of stainless steel. S.S. Safety guard rails on each side. Mounted on four castor (all with individual brakes). Provision for IV rod. Supplied completely with S.S. side railings and IV rod



Samco Medical

Website: www.samcomedical.com,