



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Stretcher Trolley

Product Code :
SAMCOMEDICAL10025A037

Description :

Samcomedical Stretcher Trolley

Technical Specification :

Stretcher Trolley Size: 1950L 650W 780H mm. Frame work made of stainless steel. S.S. Safety guard rails on each side. Basket for keeping blankets etc. Mounted on four castor. Provision for IV rod. Supplied completely with S.S. side railings and IV rod. With removable top



Samco Medical

Website: www.samcomedical.com,