



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Patient Transfer Stretcher, 2 Function

Product Code :
SAMCOMEDICAL10025A038

Description :

Patient Transfer Stretcher, 2 Function

Technical Specification :

Patient Transfer Stretcher, 2 Function Size: 1940 x 680 x 610-920mm Height adjustable by crank mechanism (920 mm-610 mm) Parallel Mattress surface movement by crank to transfer patient to Bed or table. Horizontal movement range : 400mm Aluminum $\frac{1}{2}$ Alloy Guard rails, can be fixed up and down. 150mm Central $\frac{1}{2}$ locking castors, Equipped with separate directional wheel for easy maneuvering



Samco Medical

Website: www.samcomedical.com,