



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Transfer Stretcher Cum Medical Trolley, Stainless Steel,
Adjustable Height

Product Code :

SAMCOMEDICAL10025A039

Description :

Transfer Stretcher Cum Medical Trolley, Stainless Steel, Adjustable Height

Technical Specification :

Transfer Stretcher Cum Medical Trolley, Stainless Steel, Adjustable Height Stainless Steel construction, for long life. Size: 3500 x 650 x 750 / 950 mm. Adjustable height from 750 - 950 mm by unobtrusive crank mechanism. Adjustable back section. The two stretcher base can move up / down independently and separately, and can be brought to the same level for smooth slide. Both base are locked securely before the bed slides. All wheels with individual brakes



Samco Medical

Website: www.samcomedical.com,