



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**

Samcomedical Transfer Stretcher Cum Medical Trolley

**Product Code :**

SAMCOMEDICAL10025A040

**Description :**

Samcomedical Transfer Stretcher Cum Medical Trolley

**Technical Specification :**

Samco Transfer Stretcher Cum Medical Trolley Made of high quality Aluminum alloy, PE coated, die pressed using automatic power press machines. Size: 3650 640 640 / 970 mm. Adjustable back section. Adjustable height from 640 - 970 mm. Super strong aluminum alloy tracks, increase the durability, steadiness and reliability of the moving stretcher. The two single-stretchers can move up / down independently, and can be conveniently brought to the same level for smooth slide. Both base are locked securely and automatically before the bed slides. Both single-stretchers are equipped with separate directional wheel (5th Wheel) for easy operation only by one person 125mm dia noiseless castors with simultaneous braking system which locks / unlocks 2 castors with single pedal press. International quality structure and functionality to satisfy the needs of the best hospital and Operation theaters



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),