



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Stretcher On Trolley

**Product Code :**  
SAMCOMEDICAL10025A041

**Description :**

Samcomedical Stretcher On Trolley

**Technical Specification :**

Stretcher On Trolley Size : 210L x 56W x 82H cms. Frame work made of ERW tubular pipe of 31 mm. dia mounted on four 15 cms. dia. castors (2 with brakes). Removable stretcher top made of S.S. Sheet pressed on tubular frame. Brackets on alternate ends for IV drip rod. Knockdown construction. With Oxygen Cylinder Holder



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),