



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Patient Trolley

**Product Code :**  
SAMCOMEDICAL10025A042

**Description :**

Patient Trolley

**Technical Specification :**

Patient Trolley Size : 180L x 55W x 80H cms. Frame work made of ERW tubular & rectangular pipes mounted on four 12.5 cms. dia. castors (2 with brakes). Brackets on alternate ends for IV drip rod. Supplied completely with side railings and I.V. stand. Made of chrome plated steel. Fixed Top Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),