



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Patient Trolley Removable Top (Superior)

**Product Code :**  
SAMCOMEDICAL10025A043

**Description :**

Patient Trolley Removable Top (Superior)

**Technical Specification :**

Patient Trolley Removable Top (Superior) Size : 180L x 55W x 80H cms. Frame work made of Stainless Steel Pipes mounted on four castors (2 with brakes). Brackets for IV drip rod. Supplied completely with side railings and I.V. stand. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),