



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Trolley Stainless Steel

Product Code :
SAMCOMEDICAL10025A044

Description :

Stretcher Trolley Stainless Steel

Technical Specification :

Stretcher Trolley Stainless Steel Framework made of Stainless Steel. Size (LxWxH): 190x54x85cm. Self-weight: 33 kg, Load Bearing: 159 kg. Removable Table Top. Drop down side rails. Two Diagonally opposite castor with brakes



Samco Medical

Website: www.samcomedical.com,