



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Trolley Aluminium With Adjustable Height,
Adjustable Back rest

Product Code :
SAMCOMEDICAL10025A045

Description :

Stretcher Trolley Aluminium With Adjustable Height, Adjustable Back rest

Technical Specification :

Stretcher Trolley Aluminium With Adjustable Height, Adjustable Back rest Light Weight, easy clean, anodized aluminium alloy framework with adjustable back rest. Washable, water proof & anti-fungus fabric top. Size (LxWxH): 190x55x88cm, Size of Castors : $\varnothing 160\text{mm}$. Height adjustable with crank mechanism. Supplied with I.V. Rod / Saline Rod. Range from bed surface to ground : 50cm. Self-Weight : 36 kg, Load Bearing: 159 kg



Samco Medical

Website: www.samcomedical.com,