



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electric Bed, Two Function For Home Care

Product Code :
SAMCOMEDICAL1001A073

Description :

Electric Bed, Two Function For Home Care

Technical Specification :

Electric Bed, Two Function For Home Care Two function electric bed for home care with electrically operated back rest: 0°-70° 1/2, knee rest: 0°-40° 1/2. Size: 2140mm(L)*1080mm(W)*500mm(H) All functions controlled with T-motion Power Device (Actuator), 220 - 240 V AC, 50 Hz. (110 V on request). Equipped with battery backup. Epoxy coated mild steel frame work and 4 section perforated top. Wooden side panel Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Removable high quality beautiful wooden head panel and foot panel. 125mm dia noiseless castor, all with brakes Provision for I.V. Rod on all sides of the bed



Samco Medical

Website: www.samcomedical.com,