



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Manual Bed, Two Function For Home Cares

**Product Code :**

SAMCOMEDICAL1001A075

**Description :**

Manual Bed, Two Function For Home Cares

**Technical Specification :**

Manual Bed, Two Function For Home Care Two function manual bed for home care with crank operated back rest:0°½-70°½, knee rest:0°½-40°½ Size: 2080mm(L)\*1012mm(W)\*500mm(H) Epoxy coated mild steel frame work and 4 section perforated top. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Removable high quality beautiful wooden head panel and foot panel. 125mm dia noiseless castor all with brakes. Provision for I.V. Rod on all sides of the bed



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),