



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher For Ambulance Car

Product Code :
SAMCOMEDICAL10025A046

Description :

Stretcher For Ambulance Car

Technical Specification :

Stretcher For Ambulance Car The Bed Surface of the stretcher is made of foamed cushion with adjustable backrest which makes the patient feel more comfortable. It has low frame structure, so that it can move on the ground and can also be lifted as a stretcher. The both ends of the stretcher have opening circles which are used for transporting the patients. Salvage stretcher can be put inside the stretcher. Unfolded Size (LxWxH): 190x55x42cm, Self-weight: 25 kg, Load Bearing: 159 kg. Framework made of light weight aluminium. Maximum backrest adjustable: 65½°. Suitable for the ambulances where the wheel height is less than 29cm. Supplied with HF5300 : Fixing/Locking Device for Ambulance



Samco Medical

Website: www.samcomedical.com,