



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher Cum Chair Automatic Loading For Ambulance

Product Code :
SAMCOMEDICAL10025A047

Description :

Stretcher Cum Chair Automatic Loading For Ambulance

Technical Specification :

Stretcher Cum Chair Automatic Loading For Ambulance When loading the stretcher in ambulance, the minimum height is 62cm. Bed surface of the stretcher is made of foamed cushion, the material is waterproof, the backrest is adjustable. Mechanical folding legs are controlled by the handles on both sides. Stretcher can be locked up after being aboard through the chute. Stretcher can be arbitrarily adjusted by pneumatic spring, it can converted into a chair. Framework made of aluminium. It is mainly used in limited space transporting patients, such as lift in hospital, ambulance, city road etc. Size of castor : $\varnothing 125\text{mm}$ High position : $191\frac{1}{2}\times 53\frac{1}{2}\times 82\text{cm}$, Low position : $191\frac{1}{2}\times 53\frac{1}{2}\times 25\text{cm}$, Maximum Angle of the Back : $65\frac{1}{2}^\circ$, Self-weight : 39 Kg, Load Bearing : 181 Kg. Supplied with HF5303 : Fixing/Locking Device



Samco Medical

Website: www.samcomedical.com,