



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Stretcher Automatic Loading For Ambulance

**Product Code :**  
SAMCOMEDICAL10025A048

**Description :**

Stretcher Automatic Loading For Ambulance

**Technical Specification :**

Stretcher Automatic Loading For Ambulance Light Weight, easy clean, anodized aluminium alloy framework with adjustable Back Angle 65° ± 1/2°. It has 3 step and 2 wheels of big and small according to any height of ambulance carriage bottom board. Bed surface of stretcher made of foamed cushion with safety belts which makes the patient secure and comfortable. The stretcher has safety lock. The mechanical folding legs are controlled by the handles on both sides. Only one person can send the patient to the ambulance (with this trolley) easily. Size (LxWxH): High Position: 190x55x85cm, Low Position: 190x60x26cm. Mechanical folding legs are controlled by the handles. Only one person can send the stretcher to ambulance easily. Stretcher can be locked after being aboard in Ambulance to avoid movement. Size (LxWxH): High Position: 195x55x85cm, Low Position: 195x55x25cm. Self-weight: 38 kg, Load Bearing: 181 kg. Supplied with HF5303 : Fixing/Locking Device for Ambulance



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),