



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Bowl Stand Single

Product Code :
SAMCOMEDICAL10017A001

Description :

Hospital Bowl Stand Single

Technical Specification :

Size : 35 cms S.S basin, ERW steel tubular pipe construction mounted on heavy polypropylene base with castors. Pre-treated and epoxy powder coated. Supplied complete with S.S. basin of 35 cm. dia. Knockdown construction.



Samco Medical

Website: www.samcomedical.com,