




Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name : Stretcher Ambulance With Adjustable Back Rest & Pull Along Handles	Product Code : SAMCOMEDICAL10025A049
---	--



Description :

Stretcher Ambulance With Adjustable Back Rest & Pull Along Handles

Technical Specification :

Stretcher Ambulance With Adjustable Back Rest & Pull Along Handles Light Weight, easy clean, anodized aluminium alloy framework with back rest adjustable to maximum 65° angle. This stretcher designed three-step according to the height of ambulance carriage bottom board. When loading the stretcher up to the ambulance, the lowest height is 64cm. Bed surface of the stretcher is made of hollow plastic board, light weighted, the structure is very sturdy & durable. Stretcher with safety lock, it must be locked up when transport the patients, and opened out when get on the ambulance. Mechanical folding legs are controlled by the handles on both sides. Only one person to send the patient (with this trolley) to the ambulance. The stretcher can be locked up after being aboard through the chute. It is mainly used for hospitals, ambulance and city roads carrying patients and wounded person. The size of castor: 150mm wide castors. Size (LxWxH): High Position: 195x55x85cm, Low Position: 195x55x25cm. Self-weight: 40 kg, Load Bearing: 181 kg. Supplied with HF5303 : Fixing/Locking Device for Ambulance



Samco Medical

Website: www.samcomedical.com,