



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Fixing Devices for Ambulance Stretchers

**Product Code :**  
SAMCOMEDICAL10025A050

**Description :**

Fixing Devices for Ambulance Stretchers

**Technical Specification :**

Fixing Devices for Ambulance Stretchers HF5113 & HF5115 must be supplied with HF5300 (Fixing/Locking Device for Ambulance. HF5120, HF5122, HF5141 & HF5142 must be supplied with HF5303 (Fixing/Locking Device for Ambulance)



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),