



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Anesthesia Trolley

**Product Code :**  
SAMCOMEDICAL10015A052

**Description :**

Anesthesia Trolley

**Technical Specification :**

Anesthesia Trolley Made of mild steel epoxy coated finish. Size: 650 ½ 420 ½ 1020 mm. 8 drawers for keeping medicines & goods. One Sliding shelf. Equipped with one waste bin, Douch can (SS) & oxygen holder. ABS Push handle. 4 Swivel castor all with brakes. Equipped with lock & key



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),