



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Hospital Bowl Stand Double

**Product Code :**  
SAMCOMEDICAL10017A002

**Description :**

Hospital Bowl Stand Double

**Technical Specification :**

Size : 35 cms S.S basin, Five Legs base mounted on 5 cms dia castors. With two approx. 35 cms S.S basin.  
Pre- treated can Epoxy Powder Coated. Also available in S.S. Pipe & rings



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),