



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Anesthesia Steel Trolley

Product Code :
SAMCOMEDICAL10015A054

Description :

Anesthesia Steel Trolley

Technical Specification :

Anesthesia Trolley Made of mild steel epoxy coated finish. Size: 650 420 1020 mm. 8 drawers for keeping medicines & goods. waste bin on both side of trolley. One Sliding shelf. Douch Can, S.S. & push handle. 4 Swivel castors all with brakes



Samco Medical

Website: www.samcomedical.com,