



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Sphygmomanometer Aneroid, Standard

Product Code :
SAMCOMEDICAL6001A034

Description :

Sphygmomanometer Aneroid, Standard

Technical Specification :

Sphygmomanometer Aneroid, Standard Calibrated nylon cuff, stop pin manometer, 2-tube adult size bladder, Standard latex bulb, Special vinyl zipper case, Inflation Valve with spring, Standard end valve



Samco Medical

Website: www.samcomedical.com,