



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Sphygmomanometer Aneroid Square Shaped

Product Code :
SAMCOMEDICAL6001A037

Description :

Sphygmomanometer Aneroid Square Shaped

Technical Specification :

Sphygmomanometer Aneroid Square Shaped Sphygmomanometer Aneroid Square Shaped 14cm diameter dial face adjustable for easy reading. 2 tube adult latex bladder, Std. Latex Bulb. Inflation Valve with spring, PVC Tubing



Samco Medical

Website: www.samcomedical.com,