



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Sphygmomanometer Aneroid Round Shaped - Mobile

Product Code :

SAMCOMEDICAL6001A039

Description :

Sphygmomanometer Aneroid Round Shaped - Mobile

Technical Specification :

Sphygmomanometer Aneroid Round Shaped - Mobile Sphygmomanometer Aneroid Round Shaped - Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing



Samco Medical

Website: www.samcomedical.com,