



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ophthalmoscope

Product Code :
SAMCOMEDICAL6001A046

Description :

Ophthalmoscope

Technical Specification :

Ophthalmoscope Halogen Deluxe Ophthalmoscope 2.5V Xenon Halogen Light for true tissue color and consistent, long lasting illumination 24 color coded Lenses - 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for Clear Identification of adopter setting. Ophthalmoscope handle for 2 x C size batteries (Batteries not included). Supply : Black soft pouch or hard plastic case



Samco Medical

Website: www.samcomedical.com,