



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Digital Single Channel ECG

Product Code :
SAMCOMEDICAL6001A058

Description :

Digital Single Channel ECG

Technical Specification :

Digital Single Channel ECG High-resolution thermal printing array system Frequency response is as high as 150Hz Capable of printing continuously on channel trace and annotations including lead mark and parameter such as sensitivity, paper speed and filter operation status High resolution digital filter inhibits baseline drift without affecting ECG waveform One-touch operation Concept of floating input circuit meets IEC safety standard Rechargeable battery supporting more than 50 patient examinations Lead Standard: 12 leads



Samco Medical

Website: www.samcomedical.com,