



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Directory Reader Magnifier

Product Code :
SAMCOMEDICAL6001A067

Description :

Directory Reader Magnifier

Technical Specification :

Directory Reader Magnifier Useful for reading telephone directory, computer printout, newspaper or any fine print or matter laid out in a linear format. Features a break-resistant acrylic lens with a linear scale on it to help align and measure printed matter being observed. Magnification : 2x. Field of View : 120mm x 15mm. Working Distance : Nil (in contact with print)



Samco Medical

Website: www.samcomedical.com,