



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher on Trolley

Product Code :
SAMCOMEDICAL10025A007

Description :

Stretcher on Trolley

Technical Specification :

Strong steel tubular frame work of trolley mounted on four 150 mm dia swiveling castors. Removable stretcher top. Made of aluminum sheet supported on tubular frame. Pre treated and Epoxy powder coated. Supplied with 35mm thick PU mattress. Size: 2030 mm (L) x 560 mm (W) x 820 mm (H) (Approx.) Stretcher on Trolley, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,