



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher on Trolley with Side Railings

Product Code :
SAMCOMEDICAL10025A008

Description :

Stretcher on Trolley with Side Railings

Technical Specification :

Welded tubular steel framework. Dished sheet steel top with propelling handles at each end. Mounted on 15cm dia. Castor wheels Upholstered. Supplied with 35mm thick PU mattress. With side railings. Available in : EPC Finished with S.S. Top, All EPC, All Stainless Steel Stretcher on Trolley with Side Railings, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers



Samco Medical

Website: www.samcomedical.com,