



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Deluxe Doctors Single Head Stethoscope

Product Code :

SAMCOMEDICAL6001A083

Description :

Deluxe Doctors Single Head Stethoscope

Technical Specification :

Deluxe Doctors Single Head Stethoscope Zinc alloy chest piece with nickel plated silver finish. Laser engraving brand is available. Comes with inner-spring binaural tube . Special design to get better effect through press the back of head. Special design of binaural suitable for ear structure , foam box packing. Minimum Order Qty. 3000 pcs. assorted



Samco Medical

Website: www.samcomedical.com,