



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Diagnostic set (Pin Contacting Fitting) 2.5V standard illumination for use with 2 batteries type

**Product Code :**

SAMCOMEDICAL6001A087

**Description :**

Diagnostic set (Pin Contacting Fitting) 2.5V standard illumination for use with 2 batteries type

**Technical Specification :**

Diagnostic set (Pin Contacting Fitting) 2.5V standard illumination for use with 2 batteries type  
Otoscope Head: Pin contacting fitting, Chrome plated. Fixed lens 4x in plastic mount may be removed for introductions of instruments. Reusable specula: (2.5, 3.5 and 4.5mm) included. Deluxe Ophthalmoscope head: Pin contact fitting. lens selection wheel permits choice of corrective lenses from 0 to -25 and 0 to +40 diopters single large circle aperture. Bent-arm illuminator: Tongue blade holder : Plastic nasal speculum: Chrome plated two laryngeal mirrors: Number 3 and 4 Battery Handle: Chrome plated. Rheostat for adjustment of illumination  
Require: 2x C size batteries, replacement bulb 2.5V  
Supply: Soft case or black plastic case



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),