



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Suspension Operation Theatre Light

**Product Code :**  
SAMCOMEDICAL200001A002



**Description :**

28" DIA Satellite Dome with 7 Reflectors & 22" DIA Satellite Dome with 4 Reflectors, Fully spring balanced improved high intensity light for all major surgicals procedures, Both the domes are fitted with halogen bulbs of 24V -70W each, Provided with heat absorbing and colour correcting glass filters. Satellite dome fitted with a removable type sterilizable cum focussing handle for feather touch movements. Supplied with low voltage transformer Light intensity :170000 Lux Reflector : Diachronic Colour Temp. : 4300K Temp. Increased at head area : 2~4C Finish : Oven Baked Ceiling Suspension Operation Theatre Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),