



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Air Compressor Nebulizer

Product Code :
SAMCOMEDICAL6001A103

Description :

Air Compressor Nebulizer

Technical Specification :

Air Compressor Nebulizer This equipment sprays liquid medication to fog panel by compressing air, and forms tiny particles, which flow into the throat through imbibing tube. It adopts low noise, oil-free-high-efficiency valve with no need for daily lubrication forms. It is an ideal product for family and medical unit's use
Technical Data
Maximum compressor pressure : ? 0.15Mpa Compressor free air flow : ? 10L/min Maximum nebulization rate : ? 0.20ml/min Sound level : ? 65db(A)



Samco Medical

Website: www.samcomedical.com,