



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Mercury Thermometer Clinical, Prismatic, Dual Scale,  
Economy model

**Product Code :**

SAMCOMEDICAL6001A109

**Description :**

Mercury Thermometer Clinical, Prismatic, Dual Scale, Economy model

**Technical Specification :**

Mercury Thermometer Clinical, Prismatic, Dual Scale, Economy model Oral, Prismatic type, Cylindrical bulb  
Lens front, Dual scale Measurement Range: 35-42°C & 94-108°F Accuracy  $\pm 0.1^\circ\text{C}$  &  $\pm 0.2^\circ\text{F}$  In  
plastic case



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),