



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Surgical Operation Light

**Product Code :**  
SAMCOMEDICAL200001A003



**Description :**

28in DIA Satellite Dome with 7 Reflectors, Fully spring balanced improved high intensity light for all major surgicals procedures, Both the domes are fitted with halogen bulbs of 24V -70W each, Easy manoeuvrability by semi circular handle fitted outside the dome, Angular radial and axial movement easily manoeuvrable, Supplied with low voltage transformer Light intensity :100000 Lux Reflector : Dichroic Colour Temp. : 4300K Temp. Increased at head area : 2~4C Finish : Oven Baked Surgical Operation Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),