



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Digital Clinical Thermometer

Product Code :

SAMCOMEDICAL6001A112

Description :

Digital Clinical Thermometer

Technical Specification :

Digital Clinical Thermometer Numeric display - easy to read With Beep Alert Mercury Free (environmental friendly) Unbreakable High Precision Measurement Range: 32-42°C & 89.6-109.4°F Display Resolution: 0.1°C & 0.1°F Easy Replaceable Battery Auto power off 8 minutes Last measurement memory In plastic case Digital thermometer can be used orally, rectally & under armpit



Samco Medical

Website: www.samcomedical.com,