



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Digital Blood Pressure Monitor

**Product Code :**  
SAMCOMEDICAL6001A118

**Description :**

Digital Blood Pressure Monitor

**Technical Specification :**

Digital Blood Pressure Monitor With just the touch of one button, the monitor auto inflates for quick, easy measurements. The large LCD display panel shows readings of blood pressure, pulse rate and optional SPO2 displayed simultaneously Display 2.8"½ digital LCD Power Supply 4 AA Size Alkaline Battery or optional AC/DC Adapter DP2057 - Digital Blood Pressure Monitor Easy Use With just the touch of one button, the monitor auto inflates for quick, easy measurements. The large LCD display panel shows readings of blood pressure, pulse rate and optional SPO2 displayed simultaneously. Display 2.8"½ digital LCD Power Supply 4 AA Size Alkaline Battery or optional AC/DC Adapter Configuration Adult cuff, software disc, user manual, USB line, 4 AA Size alkaline battery Features: Three-user data storage, maximum 99 records for each user. One touch auto inflate Large LCD display Three-user data system,99 records for each user Optional SPO2 function Communication with PC Clocking undisturbed in short time battery changing mmHg and kPa selectable Power off automatically when left idle for 5 minutes



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),