



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mobile OT Light

Product Code :

SAMCOMEDICAL200001A005

**Description :**

Intensity: 70,000 Lux Halogen Bulb: 24V - 70W Reflector: Diachronic Colour Temp.: 4300K Temp. Increase at Head Area: 2-4C Doma Dia: 22in Finish: Oven Baked Mobile OT Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Samco Medical

Website: www.samcomedical.com,