



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Otoscope (Pin Connection System)

Product Code :
SAMCOMEDICAL6001A135

Description :

Otoscope (Pin Connection System)

Technical Specification :

Otoscope (Pin Connection System) 2.5V Standard illumination. Chrome plated handle, 4 x magnifying glasses. Reusable specula 2.5, 3.5, 4.5mm include. Connection for pneumatic Otoscope. Battery handle for 2x AA size batteries. Rheostat for adjustment of illumination. Supply : Black soft pouch or hard plastic case



Samco Medical

Website: www.samcomedical.com,