



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ophthalmoscope with Otoscope

Product Code :
SAMCOMEDICAL6001A137

Description :

Ophthalmoscope with Otoscope

Technical Specification :

Ophthalmoscope with Otoscope Halogen Deluxe Ophthalmoscope 2.5V Xenon Halogen Light for true tissue color and consistent, long lasting illumination 24 color coded Lenses 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification of adopter setting. Supply with Otoscope handle for 2 x C Size Battery, 2.5v vacuum light (German Bulbs) High quality magnifying lens and 3 specula (size 2.5, 3.5 and 4.5) Supply : Black soft pouch or hard plastic case



Samco Medical

Website: www.samcomedical.com,