



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Professional Otoscope

**Product Code :**  
SAMCOMEDICAL6001A138

**Description :**

Professional Otoscope

**Technical Specification :**

Professional Otoscope Fiberoptic for cool, obstruction free light transmission Magnification:2.5x, Distal dia:2.5mm,4mm Bulb: 2.5V 1.6W halogen lamp Batteries: Two AA batteries or two "AA" (not included) NiMH rechargeable battery



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),