



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Surgical Mobile Light

Product Code :
SAMCOMEDICAL200001A006



Description :

Surgical Mobile Light Doma Dia: 20in, 90,000 Lux: Halogen Bulb: 24V - 150W Reflector: Diachronic Colour Temp.: 4300K Temp. Increase at Head Area: 2-4C Finish: Oven Baked Surgical Mobile Light Doma Dia: 14in, 50,000 Lux: Halogen Bulb: 24V - 100W Reflector: Diachronic Colour Temp.: 4300K Temp. Increase at Head Area: 2-4C Finish: Oven Baked Surgical Mobile Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Samco Medical

Website: www.samcomedical.com,