



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Slit Lamp Microscope

Product Code :
SAMCOMEDICAL6001A141

Description :

Slit Lamp Microscope

Technical Specification :

Slit Lamp Microscope Net slit lamp microscopes are used by ophthalmologists to directly examine the eyes of a patient under the magnification of a binocular microscope. Illumination Unit: Slit image rotation :0°½ to 180°½ Tilting Illumination :Upto 20°½ Filter Disc: Cobalt Filter, Green Filter, Yellow Filter, Two open apertures Slit diaphragm disc :Six aperture of 12,9,6,3 & 0.2mm and a wedge shaped diaphragm for infinitely variable slit lengths. Halogen Lamp: 12V,50 or 6V,20 Watts Standard Accessories 15x Eyepiece Dust Cover Plastic Breath shield Present Hruby Lens Central Screen bar Replacing Lamp (1 pc for illumination unit) Replacing Mirror Fuse Optional Accessories Applanation Tonometer Mechanical / Motorised Table



Samco Medical

Website: www.samcomedical.com,