



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Trial Lens Set

Product Code :
SAMCOMEDICAL6001A142

Description :

Trial Lens Set

Technical Specification :

Trial Lens Set NET superior Trial lens set used for various testing and diagnostic applications of eye ailments. They are appropriate for the examination of the refractive errors of the human eyes like myopia, hyperopic, astigmatism and presbyopia. They are also suitable for the testing of strabismus and color blindness. NET trial lens set are appropriate for any ophthalmology department as it is the basic tool for testing the human eyes. Technical Specifications: Spherical Lenses with handle on rims 32pairs +/- 0.12 to 20sph 20pairs +/- 0.12 to 06cyl 10 Prisms upto 12 07standard accessories Inbuilt cool illumination Supplied in beautiful wooden case with classy sunmica top and polished surface with rising tray



Samco Medical

Website: www.samcomedical.com,