



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Artificial Resuscitator (Ambu Type Bag), Silicone,
Autoclavable - Deluxe Quality (Child) and Guedel Airways

Product Code :

SAMCOMEDICAL3001A008

Description :

Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Child) and Guedel Airways

Technical Specification :

Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Child) and Guedel Airways
Ambu Type Bag Capacity 500ml. With translucent face mask size : 1, 2, & 3. With guedel airways size : 0,1 & 2.
100% latex free



Samco Medical

Website: www.samcomedical.com,