



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Artificial Resuscitator (Ambu Type Bag), Silicone,
Autoclavable - Deluxe Quality (Infant)

Product Code :

SAMCOMEDICAL3001A009

Description :

Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Infant)

Technical Specification :

Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Infant) Ambu Type Bag
Capacity 240ml. With translucent round face mask size:0. 100% latex free



Samco Medical

Website: www.samcomedical.com,