



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Surgical Operating Light Small

**Product Code :**  
SAMCOMEDICAL200001A008



**Description :**

Ceiling Surgical Operating Light Small Operation Reflector Lamp offers the illumination as compared to other standard models. Halogen Lights suitable for Major as well as Minor Surgery, fully spring balanced cool lights. Sterilizable handle has been provided in the middle of the dome for adjust the focus. A Shadow less light is suitable for all sizes of Operation Theatres as it gives maximum light output with minimum heat dissipation. Parabolic reflectors & color correcting filters. Dome Head: 550 mm x 2 Halogen Bulbs: 24V x 70 W x 8 nos. Light intensity at 1 mt: 1,00,000 Lux Height Adjustment: 914 mm Parabolic Reflector: 180 mm Color Temperature: 4200 K Control Unit: Transformer Ceiling Surgical Operating Light Small, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),